


## Electrification chemical processes

# Continuous flow reactors for electrochemistry

Organic electrochemistry is regarded a future technology for the environmentally friendly production of chemicals. There is currently a renaissance of electrochemistry, especially for the synthesis of organic compounds. It is driven by the search for “green” synthesis routes and by the emergence of novel synthetic strategies. The direct use of cheap sustainably generated (excess) electricity presents an opportunity for cost benefits over regular chemical reagents. Suitable electrochemical reactors are essential for the successful operation of electrosynthetic processes. They simultaneously need to be flexible and tailored for the desired reaction conditions.

Fraunhofer IMM has a long track in developing and realizing electrochemical reactors for a wide range of electrolytic processes. Electrochemical microreactors are characterized by their small electrode distances and high surface-to-volume ratios. Thin-gap electrochemical microreactors offer the following benefits:

- Due to the optimized design, current density and fluid distribution of high homogeneity are achieved.
- Reduction of the amount of required supporting electrolyte or even supporting electrolyte-free operation
- Large electrode surface area to volume ratio provides high space-time yield
- Integrated heat exchanger structures within a cell or reactor efficiently dissipate reaction heat, preventing hot spots and ensuring reliable temperature control.
- Small electrode distances induce low ohmic resistances and correspondingly low voltage drops.
- Enhanced selectivity by constant and uniform current density distribution, homogeneous flow and defined residence times.



*Undivided electrolysis cell for batch testing*

Engineer your technology!  
 Design your chemistry! Flow  
 electrochemistry as the  
 toolbox for building complex  
 molecules with electricity!

# Our flow electrochemistry reactors at a glance



## ELMIFLEX reactor

### Microchannel flow cell

A microstructured stacked plate reactor designed for electroorganic synthesis, especially Kolbe-electrolysis. The stainless steel plates are insulated with a PTFE coating and contain integrated heat exchangers for precise temperature control. Up to ten double cells can be combined in a stack resulting in an active electrode area of  $\leq 800 \text{ cm}^2$ . The middle electrode plate can be exchanged for a separator (diaphragm or ion exchange membrane) to create a divided cell setup.

- Material: PTFE coated stainless steel plates
- Anode: glassy carbon, BDD, graphite, stainless steel, Pt ( $5 \mu\text{m}$ ) on Ti, Ni, Ti
- Cathode: Stainless steel, Ni and Pt ( $5\text{--}10 \mu\text{m}$  coating)
- Microchannel dimensions: 56 channels (L W H  $100 \times 0.76 \times 0.15 \text{ mm}^3$ )
- Active volume:  $0.64 \text{ mL}$  per plate
- Electrode area:  $42.5 \text{ cm}^2$  per plate
- Flow rates:  $\leq 50 \text{ mL min}^{-1}$  per plate
- Operating temperatures:  $< 200 \text{ }^\circ\text{C}$
- Operating pressure:  $< 100 \text{ bar}$



ELMIFLEX reactor

## ShaPID reactor

### Zero-gap divided cell

This reactor is designed for the direct reduction of  $\text{CO}_2$  to formic acid. The zero-gap configuration is separated by an ion exchange membrane. The electrode plates are equipped with gas/liquid diffusion channels in a stainless steel current collector to contact gas diffusion electrodes and mesh electrodes.

- Material: PTFE coated stainless steel
- Anode: MMO on Ti
- Cathode: carbon GDL with electrocatalyst
- Electrode area: geometric area  $60.72 \text{ cm}^2$
- Fluid distribution microchannels: 35 per side
- Flow rates: Anode:  $0\text{--}10 \text{ mL min}^{-1}$ , Cathode:  $0\text{--}800 \text{ mL min}^{-1}$  humidified  $\text{CO}_2$
- Operating temperature:  $25\text{--}95 \text{ }^\circ\text{C}$
- Operating pressures:  $0\text{--}5 \text{ bar}$



Mechanical press for up to 5 cells



ShaPID reactor for formate synthesis from  $\text{CO}_2$



*Low-temperature electrochemical reactor*

### ELMICRYO

#### *Microchannel flow cell*

This reactor is designed for low temperature electrolysis of reactive intermediates. It uses a diaphragm or ion exchange membrane for separation of the anolyte and catholyte compartments. Temperature control is realized by contacting the reactor to a suitable cooler.

- Material: stainless steel plates coated with Pt
- Anode: Pt (5  $\mu\text{m}$  coating on stainless steel)
- Cathode: Pt (5  $\mu\text{m}$  coating on stainless steel)
- Electrode area: 12.4  $\text{cm}^2$
- Microchannel dimensions: 22 channels (L W H 44  $\times$  0.8  $\times$  0.2  $\text{mm}^3$ ), distribution and collection channels ( $\varnothing$  3 mm, L 54 mm)
- Active volume: 0.155 mL per side
- Flow rates 0.8  $\text{mL min}^{-1}$  per side
- Operating temperatures: – 80 to 60  $^{\circ}\text{C}$

### ELMISCREEN reactor

#### *Flexible divided cell*

The ELMISCREEN reactors are designed for the screening of electrocatalysts. Each half cell contains a titanium current collector equipped with gas channels making them suitable for the use of GDLs. The electrode is positioned onto the current collector. Fluid distribution of anolyte and catholyte are included in the housing. A separator (diaphragm or ion exchange membrane) can be positioned between the compartments. Different configurations include finite-gap, a zero-gap, and a high-pressure version with variable electrode areas.

- Material: insulating PEEK housing
- Anode: variable
- Cathode: variable
- Electrode area: 5–20  $\text{cm}^2$
- Cell Dimensions: 12.2  $\times$  10  $\times$  6  $\text{cm}^3$
- Flow rates: 0–50  $\text{mL min}^{-1}$  per half cell
- Operating temperatures: 20–90  $^{\circ}\text{C}$
- Operating pressures: up to 5 bar for the finite-gap and up to 80 bar for the zero-gap



*ELMISCREEN reactor for catalyst screening*



### AVATOR reactor

#### *Sandwich type parallel plate cell*

A parallel plate reactor with integrated heat exchanger suitable for synthesis or peroxod carbonate generation. High current densities are possible with sufficient cooling (up to  $1 \text{ A cm}^{-2}$ ).

- Material: insulating PEEK housing, stainless steel
- Anode: BDD, glassy carbon, Pt, other materials possible
- Cathode: stainless steel, coating with Ni, Pt possible
- Electrode area:  $2 \times 25 \text{ cm}^2$
- Typical electrode gap: 0.9 mm
- Active volume: 2.3 mL
- Flow rates: typical flow rate  $1.5 \text{ L min}^{-1}$
- Operating temperatures:  $5\text{--}80 \text{ }^{\circ}\text{C}$
- Operating pressures: 0–3 bar



AVATOR reactor



HAVANA reactor

### HAVANA reactor

#### *Sandwich-type parallel plate cell*

This reactor is designed for electroorganic synthesis consisting of stainless steel cathode plates with integrated flow distribution and heat exchangers and a middle anode plate.

- Material: Stainless steel plates
- Anode: glassy carbon, other materials possible
- Cathode: Stainless steel
- Electrode area:  $2 \times 100 \text{ cm}^2$
- Electrode gap: 1 mm
- Active volume: 10 mL
- Flow rates: optimized for  $1\text{--}6 \text{ L min}^{-1}$
- Operating temperature:  $5\text{--}120 \text{ }^{\circ}\text{C}$
- Operating pressure: tested to 3 bar

### ELMI1


#### *Undivided sandwich-type micro cell*

A compact electrochemical microreactor with integrated heat exchanger in a double cell configuration. Very small electrode gaps are possible enabling low supporting electrolyte concentrations.


- Material: PTFE-coated stainless steel, FFKM gaskets
- Anode: glassy carbon, platinum, Pb/PbO<sub>2</sub>
- Cathode: stainless steel, nickel, titanium, platinum coating, lead
- Electrode area
- Microchannel dimensions: 27 channels per plate (L W H  $29.5 \times 0.8 \times 0.1 \text{ mm}^3$ )
- Typical electrode gap: 100  $\mu\text{m}$
- Active volume: 127  $\mu\text{L}$
- Flow rates: typical flow rate  $0.1\text{--}10 \text{ mL min}^{-1}$
- Operating temperatures: up to  $200 \text{ }^{\circ}\text{C}$
- Operating pressures: up to 10 bar



ELMI1 – Electrochemical Microreactor



ECSMR Electrochemical Segmented Microreactor



### ECSMR-Screen

#### Segmented screening reactor

This cell has a segmented anode allowing for variation of the electrode area during electrolysis.

- Material: PEEK inlay, PTFE coated stainless steel housing, FFKM seals
- Anode: glassy carbon chips
- Cathode: stainless steel, coating with Ni, Pt possible
- Electrode area:  $10 \times 1 \text{ cm}^2$
- Typical electrode gap: 0.1 mm
- Active volume: 2.3 mL
- Flow rates: typical flow rate  $0.1\text{--}10 \text{ mL min}^{-1}$
- Operating temperatures: up to  $80^\circ\text{C}$
- Operating pressures: 0–3 bar

### Services and benefits

#### Our services in the field of electrochemistry in flow

We offer you a one-stop service including the entire development chain from batch to flow process testing with reactor and plant development for your dedicated flow electrochemistry equipment. Just a few examples from our application portfolio for electrochemistry in flow:

- Hofmann rearrangement [1]
- Kolbe electrolysis [2]
- Peroxod carbonate production [3]
- $\text{CO}_2$  conversion [4]
- Cation flow applications [5]
- Water electrolysis [6]

#### Feasibility studies to meet your needs

- Transfer of your batch process to continuous flow mode
- Development of your unique flow reactor for your needs
- Dedicated plant development as blueprint for your future process.

#### Contact

Dr. Thomas H. Rehm  
Division Chemistry  
Phone +49 6131 990-195  
thomas.rehm@  
imm.fraunhofer.de

Fraunhofer Institute for  
Microengineering and  
Microsystems IMM  
Carl-Zeiss-Strasse 18–20  
55129 Mainz | Germany  
www.imm.fraunhofer.de

All flyers of the  
division Chemistry  
<https://s.fhg.de/flyers-chemistry>



© Fraunhofer IMM, Mainz 2025